


Setup and Procedure


The background of the slide is a solid blue color. Overlaid on this background are several wavy, horizontal lines composed of small, dark blue dots. These lines create a sense of motion and depth, with some lines appearing more prominent than others. The lines start from the left edge and extend towards the right, with some curving upwards and others downwards.

Photodiodes

- Receives light, creates a voltage
 - Photovoltaic effect
- Must be powered
- Measures intensity of light, not wavelength



Assembling a Spectrometer



Initial design with op amp

General Design


- Powered by 9V battery
- Photodiode used to receive light
- Record intensity received by photodiode

Initial Design

- Included op amp to lower environmental interference



Final Design

- Included high-resistance resistor
- Same principle as op amp



Final circuit design

Use of an Airtrack



- Light source placed in front of airtrack
- Photodiode centered on source's center maximum
- Voltage measured with displacement from center maximum
- Higher voltage = higher intensity

Atomic Emission Setup

Blackbody Experiment (Not pictured)

- Light bulb as light source
- Diffraction gradient placed in front to separate light waves
- Data recorded with multimeter





Atomic Emission Experiment

- Similar to blackbody experiment
- Used mercury gas tube
- Start recording data before visible light

Results




Mercury Spectrum





Mercury Spectrum



- Data appears to be offset
 - Our measured peaks are longer than the known peaks
- Beam from the mercury lamp was not centered on the lens
 - Determined an angular offset of about 0.01 radians
 - Subtracted offset from each angle measurement in our data



Mercury Spectrum



Mercury Spectrum



Next Steps

The background features a solid blue color. Overlaid on this are several wavy, horizontal lines composed of small, dark blue dots. These lines create a sense of motion and depth, with some lines appearing more prominent than others, creating a layered effect.

Eliminating Error

- Isolating the spectrometer to reduce noise
- Using a better aligned lens
- Reducing oscillations in the circuit

Next Steps

- Use a rotating mirror setup
- Making it more compact
- Using software to analyze/match a spectrum